

AutoSound 9.6™ Combination Device – Automatic Hands-Free Ultrasound and Stimulation

Revolutionary patented-design allows the medical practitioner to apply ultrasound with a stationary applicator that automatically moves the ultrasound output for a more efficient and therapeutic ultrasound treatment. Combine it with four independent channels of stimulation for Hands-Free combination or concurrent ultrasound and stimulation treatments.

- Patented AutoSound design actually moves the ultrasound output for you at the ideal rate of 4cm/second
- Ultrasound output never stops moving nor moves too fast
- Focuses ultrasound treatment directly over the area – no problems with trying to treat too large an area
- Use in conjunction with the unique AutoGel Pads, made of solid conductive gel, for a perfect fit into the Gel Well to make treatment much easier, cleaner, faster and very cost-effective
- Three different size straps come standard with the device to allow the easiest strapping of the applicator to the patient
- The AutoSound is extremely safe and includes monitoring circuits plus a Patient Control Switch for added safety and peace of mind
- Comes standard with a 2cm² and 5cm² Therapy Hammer manual transducer for additional ultrasound treatment options
- Large, LCD display with touchscreen ensures ease-of-use and operation

AutoSound Ultrasound Specifications

- The AutoSound Hands-Free Applicator is capable of 1MHz, 3MHz, or 1&3 MHz Chain that alternates frequency to treat the entire tissue depth
- Capable of 2.0 W/cm² with 1MHz and 1.5 W/cm² with 3MHz or the 1&3MHz Chain
- Continuous and Pulsed Outputs (100%, 50%, 25% and 5%)
- AutoSound ERA: 3.5cm², AutoSound applicator treatment area: 14cm²
- Timer up to 20 minutes
- Nine programmable ultrasound presets
- Therapy Hammer manual transducer (2cm² and 5cm² on single soundhead)
- BNR – Maximum of 5.5:1

Stimulation Specifications

- Stimulation timer up to 99 minutes
- Four independent channels of stimulation
- None programmable stimulation presets
- Five stimulation therapy options available:
 - Quadppolar
 - Bipolar
 - Russian
 - Hi Volt (Monophasic)
 - Microcurrent

Combination Specifications

- The AutoSound Combo can provide Hands-Free Combination therapy by using the AutoSound applicator to provide both ultrasound and stimulation
- Stimulation presets are available for Combination therapy

Concurrent Stimulation and Ultrasound Treatments

The AutoSound combo allows you to provide both a Hands-Free Ultrasound and a separate stimulation treatment at the same time

Dimensions

17"W x 14"D x 7"H

Weight

19.8lbs.

Available In

120 VAC, 60Hz or 220 VAC, 50Hz