



1-888-549-4945 fax:(423)-648-7735

4120 South Creek Road Chattanooga, TN 37406

www.naimco.com www.richmarweb.com

Therasound 3.2 – Digital, Easy to use 1 MHz Ultrasound

The Therasound 3.2 is our 10cm² version of the Therasound 3 Series clinical ultrasound therapy units.

- The Therasound 3.2 has the larger 10cm² size soundhead that works very well for larger body parts.
- The 1 MHz output of the Therasound 3.2 can be used to treat the full range of applications and depths by adjusting the intensity level via the “Easy Touch” membrane pad controls.
- This ultra affordable, advanced and portable ultrasound offers both continuous and pulsed outputs.
- Additionally the 3.2 allows you to adjust Intensity, output readout(Watts/cm² or Watts) or duty cycle at any time during the treatment.
- The Therasound 3.2 also has a convenient pause feature.
- The clear, bright digital readouts and simple panel make the Therasound 3.2 easy to use.
- The 3.2’s lightweight design will work equally well in the clinic or on the road. The Therasound 3.2 comes with a standard 3 year warranty which includes the cable and transducer (most other manufacturers will only provide one year of warranty on the cable and transducer because that is the part most susceptible to breaking).

Optional Warranty Upgrade

Upgrade to a 5 Year Warranty that also includes a 2nd Day Air Loaner FREE for the entire warranty period.

Therasound 3.2 Specs

Output: 0-2.0 W/cm² (0-20 Watts)

Power Input: 110 VAC / 60 Hz, (220 VAC / 50 Hz available)

Transducers: 10cm² Effective Radiating Area (ERA)

Frequency: 1 MHz

Beam Non-uniformity Ratio (BNR): 5.5 : 1 maximum

Dimensions: 7.5” W x 7.5” D x 4”H

Weight: 5.6 lbs (2.54 kg)

Duty Cycles: 100%, 50%, 20%, 10%

Time: 0-30 minutes

Warranty: 3 Years on entire unit, including cable and transducer

Fuse: 1 amp

Line Leakage: Less than 50 mA

Therasound 3.2P Dimensions

14”W x 10.5”D x 6.5”H – 12 lbs Hard plastic case that has the ultrasound unit built into it.