

# Theratoch 4.7 Four Channel Stimulator

Touch Screen and Single "Director Dial" Controlled 4 Channel Stimulator

- Touch Screen and Single "Director Dial" Controlled 4 Channel Stimulator.
- Theratoch 4.7 has nameable, programmable QuikSets for customized treatments that can be entered and even named to meet the facility's needs allowing for quick and easy treatment setup by the whole staff
- **Stimulator** has four totally independent stimulation channels that can actually output four completely different waveforms and or treatments all at the same time.
- Upgradeable microprocessor controlled 6 Stimulation Waveforms with the ability to download any waveform to the Rich-Mar Theramini 1 satellite stimulator.
- Waveforms available on *all four channels* include:
  - **Quadpolar Interferential** (5,000 Hz carrier frequency with Cont & Chain Modes)
  - **Premodulated Bipolar Interferential** (5000 Hz carrier frequency with Cont, Surge, Alt and Chain Modes)
  - **Monophasic** (High Volt Pulsed Current twin square wave with interphase interval adjustment, with Cont, Surge, Alt and Chain Modes)
  - **Biphasic** (Symmetric square wave with phase duration and interphase interval adjustments, with Cont, Surge, Alt and Chain Modes)
  - **Russian** (2500 Hz carrier frequency with Cont, Surge and Alt Modes)
  - **Microcurrent** (Positive, negative and biphasic outputs, .3 to 1000 Hz pulse rate with Cont and Chain Modes)
- Patient safety features include low flow monitor to alarm for poor stimulation conduction and patient isolating safety loads to protect the patient from electrical shock due to power surges or equipment malfunction.
- Enclosed in an all metal cabinet
- Comes with 4 each patient lead cords and one package of electrodes
- One year warranty
- 19" W x 10"D x 5" H – 14 lbs. Available in 120 VAC, 60Hz or 220 VAC, 50 Hz

## Optional Accessories

- **Theratoch Accessory Kit**
- 8 each 4" round, reuseable carbon electrodes
- 8 each sponges for electrodes
- 8 each elastic velcro straps for electrodes

## Two Year Premium Extended Warranty

Provides two year extension of warranty that covers all parts and labor plus will provide a Loaner unit via 2<sup>nd</sup> Day Air if unit breaks. Will also arrange for pick-up and transport of broken unit back to factory for factory repair. Will then ship repaired unit back to customer and arrange for pick up and transport of loaner unit. No cost for labor, parts, loaner or shipping for the entire warranty period. Also, includes FREE upgrades of unit in the warranty period. Premium Extended Warranty is also transferable.