

Rich-Mar Winner CM2 Specifications

Dial and Button Controlled 2 Channel Stimulator and Ultrasound Combination Unit.

- **Ultrasound** has 2cm² and 5cm² soundheads on single, 4-in-1 "Therapy Hammer" transducer. Both 2cm² and 5cm² soundheads can output both 1 & 3 MHz frequencies. Can use the ultrasound while stimulation treatments are running.
- Maximum BNR of 5.5:1 with guaranteed BNR < 3.5:1 available.
- Control dial and easy touch buttons with LED indicators easily control the 4 Channel Stimulator and Ultrasound Combination Unit for quick and simple use.
- Winner CM2 has up to 25 nameable, programmable Therapy Profiles for customized treatments that can be entered and even named using the Therapy Profiler to meet the facility's needs allowing for quick and easy treatment setup by the whole staff.
- Therapy Profiles can be selected to run in either Easy Profile mode or Advanced Profile mode.
- **Stimulator** has two totally independent stimulation channels that can actually output four completely different waveforms and or treatments all at the same time. Can even use ultrasound separately while all four stimulation channels are running.
- Pulse rate options include High scan, Low scan, High-Low scan and Fixed rates.
- Can adjust pulse rates and many other parameters during the treatment.
- Microprocessor controlled 5 Stimulation Waveforms
- Waveforms available on *both channels* include:
 - **Quadpolar IFC** (5,000 Hz carrier frequency with Normal Modes)
 - **Pre-Mod IFC** (5000 Hz carrier frequency with Normal, Surge, Co-Cont, and Alt Modes)
 - **Hi-Volt** (High Volt Pulsed Current twin square wave, with Normal, Surge, Co-Cont, and Alt Modes)
 - **Russian** (2500 Hz carrier frequency with Normal, Surge, Co-Cont, and Alt Modes)
 - **Microcurrent** (.3 to 1000 Hz pulse rate in Normal Mode)
- Quad IFC has balancing feature allowing for treatment sensation focusing.
- IFC treatments have Vector massaging feature that can be turned on or off.
- Six different timing on/off modes available: 10 sec/50 sec, 10/30, 10/20, 10/10, 5/10, 5/5
- Unique channel buttons allowing for easy identification of Active treatments, real-time on or off time output and displayed parameters.
- Combination ultrasound and stimulation treatments can be used together or independently from each at the same time.
- Patient safety features include patient isolating safety loads to protect the patient from electrical shock due to power surges or equipment malfunction.
- Enclosed in an all metal clinical cabinet with stylish plastic panel frame for attractive design
- Comes with 4 each patient lead cords, one package of electrodes and 16 oz bottle of Rich-Mar Aloe Sound Ultrasound Lotion
- One year warranty
- 16" W x 9"D x 5.5" H – 13 lbs. Available in 120 VAC, 60Hz or 220 VAC, 50 Hz

Optional Accessories

Winner2 Accessory Kit

- 4 each 4" round, reusable carbon electrodes
- 4 each sponges for electrodes
- 4 each elastic velcro straps for electrodes

Two Year Premium Extended Warranty

- Provides two-year extension of warranty that covers all parts and labor plus will provide a Loaner unit via 2nd Day Air if unit breaks. Will also arrange for pick-up and transport of broken unit back to factory for factory repair. Will then ship repaired unit back to customer and arrange for pick up and transport of loaner unit. No cost for labor, parts, loaner or shipping for the entire warranty period. Also, includes FREE upgrades of unit in the warranty period. Premium Extended Warranty is also transferable.