

INTENSITY™ 7

TENS Therapy

The InTENSity™ 7 is one of the most popular, and best valued, digital TENS units on the market and is now available without a prescription. This digital TENS unit is proven to be among the most durable and easy-to-use TENS units providing drug-free pain relief for sore or aching muscles of the lower back, shoulder, arm, leg due to strain from exercise and normal activities. The upgraded digital processor is more precise and accurate than ever before. This unit is great for users seeking a digital device that is also equipped with physical amplitude knobs. FDA cleared for over-the-counter sale.

FEATURES:

- Providing Drug-Free Pain Relief for Sore or Aching Muscles of the Lower Back, Shoulder, Arm, Leg Due to Strain From Exercise and Normal Activities
- Modes: Burst, Normal, Modulation, SD1 and SD2
- Dual Channel for Use with Two or Four Electrodes
- Pulse Width and a Pulse Frequency Are Adjustable
- System Includes Device, Lead Wires, Four Self-Adhesive Reusable Electrodes, 9-Volt Battery, Hard Plastic Carrying Case, Quick Start Guide and Instruction Manual
- Compliance Meter, 60 Records Up to 999 Hours
- Low Battery Indicator
- One Year Warranty
- FDA Cleared for Over-the-Counter Sale
- Available for Private Label – Please Ask Your Rep for More Details

ORDER DETAILS:

DI0007	InTENSity™ 7
EP2020WC2-INTM	Premium White Cloth Electrodes



TECHNICAL SPECIFICATIONS

- **Power Source:** 9-volt battery
- **Channel:** Dual, isolated between channels
- **Wave Form:** Asymmetric biphasic
- **Pulse Amplitude:** Adjustable 0-100 mA, at 500 ohm load each channel
- **Pulse Rate:** Adjustable, 2-150 Hz
- **Burst Rate:** 0-5 Hz (Burst mode)
- **Pulse Width:** Adjustable, 50-300µs 10 µs/step
- **Modes:** Burst, Normal, Modulation, Strength Duration - SD1 & SD2
- **Timer:** 5 to 60 minutes with auto shutoff or continuous